

FlashDetect™ Nano

((

Fully Automated Molecular POCT System

- Result in 20~30 minutes
- Easy to Operate
- Small Footprint
- Fully Automated

All-in-One Cartridge Sample to Result

- Fully integrated NA extraction + PCR amplification
- Includes internal control
- Contamination free
- Pre-loaded with lyophilized reagents
- Room temperature storage and transportation

Vertrieb durch / distributed by:

Mast Diagnostica GmbH
Feldstraße 20
DE-23858 Reinfeld
Tel. +49 (0)4533 2007 0
Fax +49 (0)4533 2007 68
mast@mast-diagnostica.de / www.mast-group.com




TECHNICAL SPECIFICATIONS

Model	FlashDetect Nano				
Technology	Real-time PCR				
Throughput	1				
Fluorescence Channel	FAM	HEX	ROX	CY5	
Dimensions	12.5 (L) * 22 (W) * 29.5 (H) cm				
Weight	4.4 kg				



WORKFLOW



TEST MENU

	Assay	LyocartD	LyocartE*	Pack Size
	Feline Calicivirus (FCV), Feline Herpesvirus (FHV)	ð		24 Tests
	Mycoplasma felis (Mycof), Chlamydia felis (Chlf) and Bordetella bronchiseptica (Bobfla)	4		24 Tests
	Feline coronavirus (FCoV), Feline panleukopenia virus (FPV)			24 Tests
	Feline Infectious Peritonitis virus (FIPV)		٥	24 Tests
	Canine Influenza Virus (CIA), Canine Adenovirus (CADV) and Canine Parainfluenza Virus (CPIV)	ጏ		24 Tests
	Canine Distemper Virus (CDV)	<a>\bar{\bar{\bar{\bar{\bar{\bar{\bar{		24 Tests
	Afican Swine Fever Virus (ASFV)		٥	24 Tests

^{*} LyocartE assays are compatible with Coyote Flashdetect™ Nano and Flash10.

Coyote Bioscience Co., Ltd.

